



## PRODUCT SPECIFICATION SHEET

### Rascal Carbon Cruiser

FDA Class II Medical Device

Weight capacity	300 lbs
Max Speed	3.7 mph
Turning radius	31"
Unfolded Dimensions	37x24x38"
Folded Dimensions	28x24x12.5"
Front wheel/Rear wheel size	7"/8.5"
Range Per Charge	9.3 miles
Weight(without battery)	42.4 lbs
Seat Width	17.3"
Seat Depth	15.7"
Backrest Width	16"
Backrest Depth	19.2"
Motors	two - 24V/250W
Battery	24V 12Ah (LG Cell )
Battery weight	4 lbs
Controller	Brushless 35A
Max Safe Slope	9°
Seat height(From Ground)	20"
Arm rest height(From Seat)	9"
Bottom Storage Bag	Standard
Front/Rear Suspension	Yes
Carton Dimensions Size	24.8"x13"x30"